Rapid SARS-CoV-2 Antigen Test Card (Self-Test)

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Product Information 产品信息

Performance Characteristics 产品性能

- **Sample type:** nasal swab
样本类型：前鼻腔拭子
- **Time to result:** 15-20 min
结果读取时间：15-20 分钟
- **Sensitivity** 灵敏度: 96.77%
Specificity 特异性: 99.20%
Accuracy 准确性: 98.72%
- **Layperson study statistics** 非专业人员研究报告:
 - 99.10% of non-professionals carried out the test without requiring assistance
99.10%非专业人员不需要帮助可以自行完成检测及结果判读
 - 97.87% of different types of results were interpreted correctly
97.87%不同结果能被正确读取

Specifications 包装

1 Test / Kit

Catalog No.: 1N40C5-2

Kit Components:

- 1 Rapid SARS-CoV-2 Antigen Test Card (sealed foil pouch)
- 1 sterilized swab
- 1 extraction tube
- 1 sample extraction buffer
- IFU
- Tube holder on back of box

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
Test Procedures 操作步骤

①

Rotate the lid of sample extraction buffer bottle.

Caution: Open it away from your face and be careful not to spill any of the liquid.

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②


Squeeze all extraction buffer out of the bottle into the extraction tube.

Caution: Avoid touching the bottle against the tube.

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③


Identify the soft, fabric tip of the swab. Peel open the swab packaging and gently take out the swab.

Caution: Never touch the soft, fabric tip of the swab with your hands.

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④


Carefully insert swab into one nostril. The swab tip should be inserted no less than 2.5 cm (1 inch) from the edge of the nostril. Roll swab 3-4 times along the mucosa inside the nostril. Leave swab in the nostril for several seconds. Using the same swab, repeat this process for the other nostril. Withdraw swab from the nasal cavity.

Caution: This may feel uncomfortable. Do not insert the swab any deeper if you feel strong resistance or pain.

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
⑤

Place swab into extraction tube. Roll swab three to five (3-5) times. **Leave swab in extraction buffer for 1 minute.**

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
⑥


Pinch extraction tube with fingers and remove the solution from swab as much as possible. 用手指挤压管壁，从拭子上尽可能提取出溶液。

⑦

Install the nozzle cap onto the sample extraction tube tightly. 将提取管管盖安装到样品提取管上，并拧紧。

⑧


Bring the kit components to room temperature before testing. Open the pouch and remove the card. Place the card on a flat and level surface.

Caution: Once opened, the test card must be used immediately.

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⑨

Invert the extraction tube and add 3 drops (about 75 μ L, without air bubbles) of test specimen into the specimen well (S) by gently squeezing the extraction tube.

Caution: The formation of air bubbles in the specimen well (S) must be avoided.

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⑩



Read the results at 15-20 minutes.

Caution: Results after 20 minutes may not be accurate.

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Interpretation of Results 结果判读

 <p>Positive</p>	<p>Positive: If two colored bands appear with one colored band in the Control Zone (C) and another in the Test Zone (T) within 15-20 minutes, the test result is positive.</p> <p>Caution: No matter how faint the colored band is in the Test Zone (T), the result should be considered as positive.</p> <p>阳性: 若 15-20 分钟内质控线 (C) 和检测线 (T) 都出现紫红色条带, 表示检测结果有效且为阳性。</p> <p>注意事项: 不论检测线 (T) 显色强度有多弱, 结果都应视为阳性。</p>
 <p>Negative</p>	<p>Negative: If one colored band appears in the Control Zone (C) and no colored band appears in the Test Zone (T) within 15-20 minutes, the test result is negative.</p> <p>阴性: 若 15-20 分钟内仅在质控线 (C) 出现紫红色条带, 检测区 (T) 未出现紫红色条带, 则检测结果有效且为阴性。</p>
 <p>Invalid</p>	<p>Invalid: If no color line appears in the control area (C) within 15-20 minutes, the test is invalid. Repeat the test with a new test card.</p> <p>无效: 若 15-20 分钟内质控线 (C) 未出现紫红色条带, 表示检测无效, 应重新检验。</p>

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
ZERTIFIKAT ◆ CERTIFICATE ◆ 證書 ◆ CERTIFICADO ◆ CERTIFICAT		Benannt durch/Designated by Zentralstelle der Länder für Gesundheitsschutz bei Arzneimitteln und Medizinprodukten www.zlg.de ZLG-BS-245.10.07		
	EC Certificate			
	EC Design-Examination Certificate Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III (6) (Devices for self-testing)			
	No. V9 061317 0006 Rev. 00			
	Manufacturer:	Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road Jimei North Industrial Park 361021 Xiamen, Fujian PEOPLE'S REPUBLIC OF CHINA		
	Product:	In Vitro diagnostic devices for self testing		
	Model(s):	Rapid SARS-CoV-2 Antigen Test Card		
	Parameters:	Model Name:	Model No.:	
		Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-2	
		Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-4	
	Rapid SARS-CoV-2 Antigen Test Card	REF 1N40C5-6		
The Certification Body of TÜV SÜD Product Service GmbH declares that a design examination has been carried out on the respective devices in accordance with IVDD Annex III (6). The design of the devices conforms to the requirements of this Directive. All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:V9_061317_0006_Rev_00				
Report No.:	713210321			
Valid from:	2021-04-01			
Valid until:	2022-05-26			
Date,	2021-04-01			
				
	Christoph Dicks Head of Certification/Notified Body			
Page 1 of 1 TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123 TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany				
				

Declaration of Conformity

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Declaration of Conformity	
Manufacturer	Xiamen Boson Biotech Co., Ltd. 90-94 Tianfeng Road, Jimei North Industrial Park, Xiamen, Fujian 361021, P. R. China.
European Representative	Lotus NL B.V. Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands.
Product	Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-2 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-4 Rapid SARS-CoV-2 Antigen Test Card, REF: 1N40C5-6
Classification	Self-Test IVD
Confirmaty Assessment Route	Directive 98/79/EC on In Vitro Diagnostic Medical Devices (IVDD), Annex III section 6
We herewith declare that above mentioned products meet the provisions of the council Directive 98/79/EC for medical devices. All Supporting Documentation is Retained under the premises of the manufacturer. We have sole responsibility for issuing the Declaration of Conformity.	
Standard Applied	ISO 13485:2016 EN ISO 14971:2012 EN ISO 18113-1:2011 EN ISO 18113-4:2011 ISO 15223-1:2016 EN 62366-1:2015 EN 13612:2002 EN ISO 23640:2015
Notified Body	TÜV SÜD Product Service GmbH, Ridlerstraße 65, 80339 Munich, Germany. Identification number: 0123
(EC) Certificate(s)	No. V9 061317 0006 Rev. 00
Start of CE-Marking	2021-04-01
Place, Date of Issue	Xiamen, 2021-04-01
Signature	 (Signed By Boson Representative) Name: Changgong Zhang Title: General Manager
	

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Authorization Acknowledgement from China Chamber of Commerce for Import & Export of Medicines & Health Products 中国医药保健品进出口商会白名单

Source 网址: <http://www.cccmhpie.org.cn/Pub/6325/178081.shtml>

取得国外标准认证或注册的医疗物资企业清单					
Name List of Medical Devices and Supplies Companies with Certification/Authorization from other Countries					
序号	生产企业	统一社会信用代码	国外注册认证情况	省份	产品型号
No.	Company	Uniform Social Credit Code	Status of Certification / Authorization in Other Countries	Province	Product Model
60	厦门市波生生物技术有限公司 Xiamen Boson Biotech Co., Ltd	91350200705468594R	CE	福建	2019-nCoV IgG/IgM Combo Test Card Rapid SARS-CoV-2 Antigen Test Card Rapid SARS-CoV-2/Influenza A+B Antigen Duo Test Card Anti-SARS-CoV-2 S Protein Antibody Quantitative Test SARS-CoV-2 Neutralizing Antibody Test Card SARS-CoV-2 Neutralizing Antibody Quantitative Test SARS-CoV-2 Antigen Saliva Test Card Anti-SARS-CoV-2 S Protein Antibody Test Card

Registration Completed in Individual Countries

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